

HAAKE Viscotester iQ RheoApp Start Guide


The HAAKE Viscotester iQ RheoApp software is a tool for editing (internal) Jobs, setting up the (internal) User Management system and editing the configuration of the HAAKE Viscotester iQ user interface. The software runs directly from the USB flash drive it is supplied on.



Setting up the USB flash drive

Starting from RheoApp version 2.4.4.18, RheoApp does not need to be paired with a specific HAAKE Viscotester iQ rheometer anymore.

The HAAKE Viscotester iQ RheoApp flash drive can be used as delivered, with *any* HAAKE Viscotester iQ, a setup and a license key are *not* needed.



Working with the software and the USB flash drive

The handling of the HAAKE Viscotester iQ USB flash drive for transferring Jobs and data from and to the HAAKE Viscotester iQ rheometer, as well the operation of the HAAKE Viscotester iQ RheoApp software is described in detail in the HAAKE Viscotester iQ user guide.

Contents of delivery

Part	Number
HAAKE Viscotester iQ USB flash drive	098-5065
with HAAKE Viscotester iQ RheoApp software	
This Start Guide	006-2377