


Example Application Document

In a Job, replace this document with a SOP document, a document describing this Job, a document describing the application or any other document that is useful.


This example document describes the functionality of the Application information area of the Job-Editor.

Application information area in the Job-Editor

In the Application information area a short text describing the Job is always visible, unless the area is collapsed to its minimum height.



❖ To show the Application information

1. Click the information  button.


❖ To open the PDF document



1. Click the PDF  button.

❖ To edit the Application information


1. Click the **Edit** button.

❖ To collapse the Application information area

1. Click the collapse  button.

When no Long text is available the Information  button is not displayed. When no Document file name is specified the PDF  button is not displayed.


The application information edit dialog

In the Application information edit dialog the Short text and the Long text for the Application information can be edited. These texts are stored in the Job file. If desired a media (image and/or video) file can be selected for display together with the long text in the Application information. A Document (PDF) file can be selected for display from the Job-Editor using the PDF  button.

Note Only the path to the media and document files are stored in the Job, not the actual files.

Most of the example Jobs that are part of the RheoWin installation contain predefined Long text as well as images or videos and the PDF buttons links to this document.


Note Although both an image file name and a video file name can be selected (and stored in the Job file), either the image or the video (or no media at all) is displayed in the Application information display window.



The default directory for the media and document files is C:\Users\Public\Documents\Thermo\RheoWin\Media it is recommended to always use this directory for files that are used in Jobs.

The application information display

When both a Long text and a media (image or video) file are defined the Application information is displayed as in the following figure.



When no media file is selected the Application information is displayed in a simple message window.